

KKM / KKM Puzzle / KKM LongLine









the cosy mat classic since 2002

- the highly wear resistant honeycomb profile creates excellent and permanent lying comfort
- integrated soft slope to rear edge promotes drying off
- proven pebbled surface
- KKM LongLine: installation independent of divider arrangement



the sealing lips have proven themselves in practice



sealing lips on lower side

options:



with straight edges



with KK-Poly-T-Bar e.g. for stall equipment with mushroom dividers or especially when using very fine litter



with puzzle shaped edges



as continuous system KKM LongLine

Dimensions:

thickness: 30 mm

KKM (straight edges):

110/115/120/125/130 cm

length:

170/183 cm

KKM Puzzle: width: length:

110/120/125 cm

170/183 cm

KKM LongLine: length:

2.5 - 50 m

width:

170/183/200 cm

in 10 cm increments

Installation*: 2 - 3 fixings / lying place

Accessories:

KK-Poly-T-Bar for covering joint widths of 40 - 80 mm or for closure at side edges, length: 120/165/180 cm,

installation: 4 - 5 screws / profile

ergoBOARD / maxiBOARD - brisket boards (see page 16)

stPlease take notice of the corresponding KRAIBURG specifications and installation instructions with useful advice and tips. You will obtain these at your local KRAIBURG dealer or at www.kraiburg-elastik.com.